Instructions for use
Endonasal postsurgical activator Nose-Fit™

Destination

Nose-Fit™ is an endonasal retainer aiming to ensure nasal breathing, shape the nasal dome and maintain the septum in a median position. Nose-Fit™ is designed to prevent deformation during scarring after rhinoplasty, septoplasty and rhinoseptoplasty procedures.

Description

Nose-Fit™ was designed using 3D-modeling techniques. The shape of the activator closely matches the normal anatomy of the internal nasal cavities and assures an optimal function.

Nose-Fit™ is comfortable for the patients. No fixation device is required as it is securely anchored in the nasal cavity due to its specifically designed shape. The largest part of the activator stays inside the nose. It is thus a discrete device that can be worn without interfering with a normal lifestyle. Nose-Fit™ is made of two tubes inserted in the nasal cavities and connected by an arch-like structure. We used high quality medical silicone in order to minimize the risks of allergic reactions and nasal mucosal damage. Nose-Fit™ is available in several standard sizes. The appropriate size has to be chosen by the surgeon using the provided table. The measurement technique is explained in Figure 2.

Before using

Before use, the Nose-Fit™ retainer has to be sterilized in boiling water for 5 minutes or disinfected using chlorhexidine. Before each insertion, the retainer has to be moisturized using boiled water or sterile saline in order to protect the nasal mucosa.

Before use, the Nose-Fit™ retainer may need bespoke tailoring. This procedure is conducted by the surgeon and can consist in reducing the length of the tubes and/or modifying the arch-shaped region (Fig. 3). These modifications do not interfere with the function of the retainer.

Activator placing

The position of the activator in the nose is shown in Figure 4. The device has to be carefully inserted simultaneously in both nostrils until the whole

tubes are inside the nose and the connecting arch is at the level of the columella. Smooth forceps may be of use in the insertion process.

If Nose-Fit™ is to be used in children under 2 years of age, it is necessary to fix the retainer using medical tape to the dorsal part of the nose when the young patient cannot be looked after closely (for instance during sleep).

The activator is designed for constant use during the whole post-surgical recovery and scar healing period (1–6 months). The Nose-Fit™ activator must be replaced in case of damage or loss.

Care of activator

In order to ensure optimal hygienic conditions when using Nose-Fit™, take out the activator at least once a day and clean it up inside and outside using running water. Cleaning of the nasal cavity has to be performed according to the physicians indications when the retainer is taken off. The retainer has to be put back in functional position immediately after cleaning.

It is recommended to instillate drops of isotonic saline solution into the retainer openings several times a day in order to prevent dryout and crusting.

Trouble shooting

If inserting the retainer is difficult or if it is loose once fully inserted, the size may be wrong. In this case, get in touch with your surgeon for replacement.

Important

Please use the Nose-Fit™ activator under medical prescription only! Your doctor may define the conditions and the length of the treatment. The recovery period following surgery has to be supervised by your surgeon. The consequences of surgery and of the potential use of the activator are the responsibility of your surgeon only.

Nose-Fit™ is a registered trade mark. All rights reserved.

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